

INFECTIOUS DISEASE

AIDS

Anti-HIV TETRA ELISA (anti-HIV-1/2)
N0001/96 Wells - N0002/480 Wells

Intended Use	For the qualitative determination of HIV-1/HIV-2 Antibodies.
Label	HRP
Method	ELISA (Break-away Wells)
Controls Included	Two Controls (positive, negative)
Incubations	1 h + 1 h/37°C.
Chromogene/Substrate	TMB
Sample Types	Serum/Plasma
Sample Size	150 µL
Comments	Not available in all countries