

BONE AND MINERAL METABOLISM

BONE MARKERS - RIA

N-tact® Osteo SP
32300/100 Beads

Intended Use	For the quantitative determination of intact human Osteocalcin (1-49).
Label	125 I
Method	IRMA Beads
Controls Included	Two Levels
No. Standard	6; 0-200 ng/mL
Detection Limit	0.2 ng/mL
Incubations	2 h/RT, shaking
Feature & Benefits	Excellent dilution linearity. Small sample size. Wide dynamic range
Sample Types	Serum
Sample Size	20 µL
Radioactivity	1110 kBq - 30 µCi/100 tests