

ENDOCRINOLOGY

GASTROENTEROLOGY - RIA

GammaDab Gastrin CA1570/125 Tubes

Intended Use	For the quantitative determination of Gastrin.
Label	125 I
Method	RIA Double Antibody
Controls Included	One Level
No. Standard	6; 0-1000 pg/mL
Detection Limit	6 pg/mL.
Incubations	1 h/RT
Sample Types	Serum
Sample Size	100 µL
Radioactivity	74 kBq - 2 µCi/125 tests

TRYPSIK P2573/100 Tubes

Intended Use	For the quantitative determination of Trypsin-like Immunoreactivity (TLI).
Label	125 I
Method	RIA Double Antibody
Controls Included	One Level
No. Standard	7; 0-400 ng/mL
Detection Limit	2.5 ng/mL.
Incubations	1 h/RT
Sample Types	Serum/Plasma
Sample Size	100 µL
Radioactivity	57 kBq - 1.55 µCi/100 tests