

AUTOIMMUNITY

TISSUE SECTIONS - IFA

CT-3 FLUORO-KIT Rat Kidney, Liver, Stomach 1740/10 x 6 Wells

Intended Use	For the detection of ANA, AMA, ASMA, LKM, APCA Antibodies.
Method	IFA
Controls Included	Homogenous, AMA, ASMA Positive, Negative Controls
Incubations	30 min + 30 min/RT
Feature & Benefits	Triple Substrate. Superior Hydrophobic Mask eliminates blotting.
Conjugates	FITC
Reading	Fluorescent Microscope
Sample Types	Serum

CT-2 FLUORO-KIT Rat Kidney, Stomach 1741/10 x 6 Wells

Intended Use	For the detection of ANA, ASMA, APCA Antibodies.
Method	IFA
Controls Included	AMA, ASMA Positive, Negative Controls
Incubations	30 min + 30 min/RT
Feature & Benefits	Dual Substrate. Superior Hydrophobic Mask eliminates blotting.
Conjugates	FITC
Reading	Fluorescent Microscope
Sample Types	Serum

RK FLUORO-KIT Rat Kidney 1960/10 x 4 Wells - 1730/10 x 12 Wells 1730/10 x 8 Wells

Intended Use	For the detection of AMA Antibodies.
Method	IFA
Controls Included	Homogenous, AMA Positive, Negative Controls
Incubations	30 min + 30 min/RT
Feature & Benefits	Superior Hydrophobic Mask eliminates blotting.
Conjugates	FITC
Reading	Fluorescent Microscope
Sample Types	Serum

RS FLUORO-KIT Rat Stomach

1970/10 x 4 Wells

Intended Use	For the detection of ASMA, APCA Antibodies.
Method	IFA
Controls Included	Homogenous, AMA Positive, Negative Controls
Incubations	30 min + 30 min/RT
Feature & Benefits	Superior Hydrophobic Mask eliminates blotting.
Conjugates	FITC
Reading	Fluorescent Microscope
Sample Types	Serum

RL FLUORO-KIT Rat Liver

1930/10x4 Wells - 1890/10x8 Wells - 1770/20x8 Wells

Intended Use	For the detection of Antinuclear Antibodies.
Method	IFA
Controls Included	Homogenous Positive, End point Reference Negative Controls
Incubations	30 min + 30 min/RT
Feature & Benefits	Superior Hydrophobic Mask eliminates blotting
Conjugates	FITC
Reading	Fluorescent Microscope
Sample Types	Serum

EMA FLUOR

N0035/10 x 5 Wells - N0049/10 x 10 Wells N0049/10 x 10 Wells

Intended Use	For the detection of Antithyroid Antibodies.
Method	IFA
Controls Included	EMA Positive, Negative Controls
Incubations	30 min + 30 min/RT
Feature & Benefits	Superior Hydrophobic Mask eliminates blotting.
Conjugates	FITC
Reading	Fluorescent Microscope
Sample Types	Serum

ICA FLUOR

N0036/10 x 5 Wells-N0074/10 x 10 Wells N0074/10 x 10 Wells N003610/10 x 5 Well

Intended Use	For the detection of Islet Cell Apparatus Antibodies (ICA).
Method	IFA
Controls Included	ICA Positive, Negative Controls
Incubations	30 min + 30 min/RT
Feature & Benefits	Superior Hydrophobic Mask eliminates blotting.
Conjugates	FITC
Reading	Fluorescent Microscope
Sample Types	Serum

ASA FLUOR

N0037/10 x 5 Wells - N0073/10 x 10 Wells N0073/10 x 10 Wells

Intended Use	For the detection of Anti-Skin Antibodies (ASA).
Method	IFA
Controls Included	ASA Positive, Negative Controls
Incubations	30 min + 30 min/RT
Feature & Benefits	Superior Hydrophobic Mask eliminates blotting.
Conjugates	FITC
Reading	Fluorescent Microscope
Sample Types	Serum