

AUTOIMMUNITY

COMPLEMENT ACTIVATION - ELISA

Complement Activation EIA (CAE) Kit 6200/96 Wells

Intended Use	For the measurement of Total Classical Complement Activity.
Label	HRP
Method	EIA
Controls Included	Two Levels-High & Low
No. Standard	1
Incubations	1 hr +1 hr + 30 min/37°C
Chromogene/Substrate	TMB
Feature & Benefits	96 tests. Simple and Objective. Proven Clinical Utility.
Reading	Spectrophotometric
Sample Types	Serum
Sample Size	5 µL
Comments	In Vitro Diagnostic Use