

## SPQ PROTEIN

---

### FITC CONJUGATES

Lipoprotein (a)  
82950/2 mL

<b>Intended Use</b>	The high titer and avidity of this antiserum allows for a wide range of applications including nephelometry and qualitative and quantitative gel techniques.
<b>Method</b>	Nephelometric/Turbidimetric
<b>Feature &amp; Benefits</b>	Antisera produced in goat. Liquid reagent, stable for up to 2 1/2 years at 2-8°C. Quality controlled through analysis by rate nephelometry.
<b>Comments</b>	Analyte Specific Reagent*