

ENDOCRINOLOGY

GASTROENTEROLOGY - RIA**GammaDab™ Gastrin**
CA1570/125 Tubes

Intended Use	Quantitative determination of Gastrin levels.
Label	125 I
Method	RIA Double Antibody
Controls Included	One Level
No. Standard	6
Incubations	70 min/RT
Sample Types	Serum
Sample Size	100 µL
Comments	For In Vitro Diagnostic Use.

TRYPSIK
P2573/100 Tubes

Intended Use	Quantitative determination of Trypsin-like immunoreactivity (TLI).
Label	125 I
Method	RIA Double Antibody
Controls Included	One
No. Standard	7
Incubations	1 hr/RT
Sample Types	Serum/Plasma
Sample Size	100 µL
Comments	For In Vitro Diagnostic Use.