

ENDOCRINOLOGY

HEMATOLOGY - RIA

Vitamin B12 / Folate
06B257117/100 Tubes

Intended Use	Simultaneous quantitative determination of Vitamin B12 and Folate.
Label	125 I, 57CO
Method	RIA Dual Isotope
Controls Included	Not Included
No. Standard	6
Incubations	15 min + 10 min + 1 hr/RT
Sample Types	Serum/Plasma/Whole Blood
Sample Size	200 µL
Comments	For In Vitro Diagnostic Use.