

INFECTIOUS DISEASE

INFECTIOUS DISEASE - IFA

CMV-vue HRP
24300/8 x 8 Wells

| | |
|-----------------------------|---|
| Intended Use | For the Detection of early structural protein (pp65) of Cytomegalovirus |
| Label | HRP |
| Method | Direct Staining |
| Controls Included | Each slide containing one positive control well and one negative control well |
| Incubations | 3 hrs (sample preparation); 1.5 hrs (assay) |
| Chromogene/Substrate | TMB |
| Sample Types | Leukocyte Isolate |
| Sample Size | 3 mL |
| Comments | In Vitro Diagnostic Use |