

LIAISON®**OTHER INFECTIOUS DISEASES****LIAISON® Treponema Assay**
310480

Intended Use	For the qualitative determination of total antibodies to <i>Treponema pallidum</i> .
Material Provided	Reagent Integral for 200 determinations. Diluent included.
Method	1-step chemiluminescent immunoassay using directly coated magnetic microparticles.
Calibrators	Included
Incubations	30 min
Sample Types	Serum

LIAISON® Borrelia burgdorferi
310870

Intended Use	For the qualitative determination of IgG and IgM antibodies to <i>Borrelia burgdorferi</i>
Material Provided	Reagent integral for 100 determinations. Diluent included.
Method	2-step chemiluminescent immunoassay using directly coated magnetic microparticles.
Calibrators	Included
Incubations	20 min
Sample Types	Serum

LIAISON® VZV IgG
310490

Intended Use	For the qualitative determinations of IgG antibodies to Varicella-Zoster virus.
Material Provided	Reagent Integral for 50 determinations. Diluent included.
Method	2-step chemiluminescent immunoassay using directly coated magnetic microparticles.
Calibrators	Included
Incubations	20 min
Sample Types	Serum