

## AUTOIMMUNITY

---

### NDNA - IFA

nDNA FLUORO-KIT *Crithidia luciliae*  
1790/10 x 5 Wells-1860/10 x 10 Wells 1860/10 x 10 Wells

<b>Intended Use</b>	For the detection and titration of circulating Native DNA (nDNA) Antibodies.
<b>Method</b>	IFA
<b>Controls Included</b>	nDNA Positive, Negative Controls
<b>Incubations</b>	30 min + 30 min/RT
<b>Feature &amp; Benefits</b>	Superior Hydrophobic Mask eliminates blotting.
<b>Conjugates</b>	FITC
<b>Reading</b>	Fluorescent Microscope
<b>Sample Types</b>	Serum