


Confidence in Your Results

LIAISON® Rubella IgG (code 310720)

LIAISON® Rubella IgM (code 310730)

A unique selection of raw materials

- * Raw materials produced under controlled conditions
- * DiaSorin uses a sucrose-gradient purified viral particle to assure highest specificity



Assay format ensures reliable results

- * Superior sensitivity allows detection of early infection
- * Superior diagnostic and analytical specificity allows curbing of interference from rheumatoid factor (RF)
- * Luminescence technology in combination with magnetic microparticles
- * Wide measuring ranges
- * Quantitative determination of IgG (IU/mL) and IgM (AU/mL)
- * Controlled reagent cooling and incubation conditions

Ease of use

- * Full automation makes your daily routine convenient and easy
- * Barcoded samples and reagents
- * Two-point recalibration, stable for 4 weeks
- * Flexible assay protocols
- * Continuous reagent inventory
- * Calibrators included
- * Ready-to-use reagent cartridge

Flexibility enables quick results

- * High throughput:
 - * Rubella IgG: 90 results/hour
 - * Rubella IgM: 50 results/hour
- * Time to first result:
 - * Rubella IgG: 35 min
 - * Rubella IgM: 45 min
- * Stored master curve
- * Small sample volume



DiaSorin S.p.A.
Via Crescentino
13040 Saluggia (VC) - Italy
Tel. +39.0161.487093 - Fax: +39.0161.487628

E-mail: info@diasorin.it
www.diasorin.com

DiaSorin