


HR

LIAISON® mureX Anti-HEV DIAGNOSTIC SOLUTION

The combined testing solution
for the detection of antibodies
to Hepatitis E virus

HEPATITIS AND RETROVIRUSES



FOR OUTSIDE US
AND CANADA ONLY

DiaSorin

The Diagnostic Specialist

LIAISON® MUREX Anti-HEV IgG & IgM

The first fully automated chemiluminescence immunoassays for the diagnosis of HEV infection and the selection of organ, tissue and cells post-mortem donors.

Clinical overview

- Hepatitis E virus (HEV) is the most common cause of acute viral hepatitis worldwide¹
- HEV infection is mostly asymptomatic and self-limiting but it can lead to fulminant hepatitis in pregnant women and patients with pre-existing liver disease
- The concentration of anti-HEV IgG antibodies could be useful in determining which level of IgG would be protective after natural infection or vaccine administration²
- According to EASL guidelines, all patients with hepatitis should be tested for HEV, as part of the first-line virological investigation³
- Most recent evidence on HEV transmission after transplantation suggests the need for a systematic screening of donors, particularly in endemic regions⁴

The value of the DiaSorin solution

Boost the global response to HEV with DiaSorin solution.

- Provide reliable results choosing the good performances of the LIAISON® MUREX HEV assays
- Enhance the diagnostic value of your results providing the quantification of IgG anti-HEV antibodies
- Be recognized as a reference Center for the organ donors selection by using the widest fully automated, chemiluminescent panel validated for the use on samples collected from dead donors
- Improve your capability to support the clinical patients management thank to the differentiation of five different viral hepatitis infections
- Improving the efficiency and the workflow of your laboratory through the consolidation on the LIAISON® XL platform of the widest automated menu for the diagnosis of Infectious Diseases and Viral Hepatitis

The assays features

Clinical Performances:

Anti-HEV IgM: 97.93% Specificity; 98.44% (PCR positive) Sensitivity
Anti-HEV IgG: 99.22% Specificity; 97.92% Sensitivity

Expression of results:

Anti-HEV IgG: quantitative (IU/mL)
Anti-HEV IgM: qualitative (index)

Throughput:

LIAISON® XL: 171 test/h
LIAISON®, LIAISON® XS: 86 test/h

Time to first result:

35 minutes

Sample Type:

Human Serum
Human Plasma (Li-heparin, Na-heparin, K2-EDTA, K-oxalate, ACD, CPDA)

Reagent, Controls and Calibrators:

Ready to use

Ordering Information

Description	Kit configuration	Part number
LIAISON® Murex Anti-HEV IgG	100 tests	311270
LIAISON® Murex Control Anti-HEV IgG	2 levels/2 vial each	311271
LIAISON® Murex Anti-HEV IgM	100 tests	311280
LIAISON® Murex Control Anti-HEV IgM	2 levels/2 vial each	311281

References

1. Mehnaaz S Khuroo et al., J Viral Hepat. 2016 Feb;23(2):68-79
2. Kamar N. et al., Clin Microbiol Rev 2014;27:116-138
3. EASL Clinical Practice Guidelines on hepatitis E virus infection. 2018
4. Pourbaix A. et al., Transpl Infect Dis. 2017 Feb;19(1)

Available on LIAISON® systems

Product availability subject to required regulatory approval

Diasorin S.p.A

Via Crescentino SNC, 13040 Saluggia (VC) - Italy
www.diasorin.com